GerauchsinformationDE

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DESCRIPTION
OArt+ is a sodium hyaluronate gel dosed to 1 %, 1,6 %, 2,2 % or 2,5 % which is obtained by bacterial fermentation. It is as a viscoelastic gel to supplement the intraarticular connective tissue. It is sterile and apyrogen.

COMPOSITION
1 ml of the product contains 10 mg, 16 mg, 22 mg or 25 mg of sodium hyaluronate with fermentative origin, low phosphate tampon (sterile, isotonic, apyrogen, pH 7.2), and water for injection.

INDICATIONS
OArt+ is injected into the joint to relieve pain and improve articular mobility associated with degenerative and trauma knee lesions and other synovial joints.

CONTRA-INDICATIONS
The product should not be used for patients
- with a known hypersensitivity reaction to any of the components of the product
- with septic arthritis
- with infections of the skin or dermatological diseases at the injection side
- given coagulation inhibitors such as Macumar

APPLICATION WITH CAUTION
Since septic arthritis is a serious side effect, please observe all usual precautions as for surgical interventions.
The product should be injected into the articulation. Intra-vascular or intra-tissular injection should be avoided.
Since no clinical data are available on the use of hyaluronic acid in children or pregnant woman, these patients should not be treated.
If packaging is damaged product should not be used.
The product has to be used before the expiry date indicated on the packaging.

SIDE EFFECTS
No sign of bioincompatibility has been detected throughout biocompatibility testing of the product.
Some local transient reactions may occur after intra-articular injection, such as redness, swelling, heat and pain. The local application of ice for 5-10 minutes may reduce these effects.

INTERACTIONS WITH OTHER MEDICAL PRODUCTS
No known incompatibility with other intra-articular injectable products has been reported up to now.
Simultaneous intake of analgesics or anti-inflammatory drugs by oral route may be useful during the first days after treatment.
The product must not be used in association with quaternary ammonium compounds.
The product must not be resterilised or reuse because this could lead to a risk of infection and/or cross-infection and/or malfunction of the product.

ADMINISTRATION AND DOSAGE
The syringe is prefilled. The Luer lock closure is removed and a suitable sharp edged tip sterile needle is screwed on the Luer lock. Use of a single use needle is recommended.
Remove air bubbles from the syringe before injection.
The skin is disinfected with local antiseptic.
The product is injected in the joint, once a week, 1-5 times.
The beneficial effects of 1-5 sequential treatments are maintained during at least 6 months.
In the case of intra Articular swelling, fluid should be removed by aspiration before the hyaluronic acid injection.
Depending on the size of the joint: 1, 3 or 5 injections of 1ml or 2ml may be injected intra-articular.
- The prefilled syringe of 1ml is the volume adapted for small joints.
- The prefilled syringe of 2ml is the volume adapted for the knee, hip, shoulder and ankle.
- The prefilled syringe of 4.8 ml is the volume adapted for the knee only.
The recommended treatment regimen with the mono dose of 4.8 ml for the knee is one injection, which may be repeated 6 months later based on the clinical experience and the pain relief obtained by the patient.

PROPERTIES AND MODE OF ACTION
The synovial fluid contains physiologically hyaluronic acid which is responsible for viscoelastic properties allowing, through lubricant properties, painless joint movements. It also participates to the nutrition of the cartilage. Intra-articular supplementation of the synovial fluid by hyaluronic acid injections improves the viscoelastic properties of the synovial fluid. After 1 - 5 intra-articular injections, joint mobility is improved, and the pain associated with degenerative lesions is relieved.

PACKAGING
a) 1,0 ml of intra articular gel contained 10 mg of sodium hyaluronate in prefilled syringe. 1 syringe are placed in sterile blister packaging. 1 blister together with instruction for use in carton pack.
b) 2,0 ml of intra articular gel contained 44 mg of sodium hyaluronate in prefilled syringe. 1 syringe are placed in sterile blister packaging. 1 blister together with instruction for use in carton pack.
c) 4,8 ml of intra articular gel contained 120 mg sodium hyaluronate in prefilled syringe. 1 syringe are placed in sterile blister packaging. 1 blister together with instruction for use in carton pack.

Exclusively use by medical doctors.

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